Innovative Powder Strength Tester



Sample Size	Consolidation Pressure	Quick Results
The old way:	The old way:	The old way:
50 to 600 grams per test	1000 Pa (1kPa) or more	as many as 5 tests to get one point
The new way:	The new way:	The new way:
0.06 to 0.1 gram per test	from 10 Pa to 50 kPa	one test = one point



Comes Complete With

- pre-programmed laptop
- bonus testing cell
- instruction manual with demonstration video
- one-year limited warranty on parts and labor
- micro-spatula for easy fill

Manufactured by:



Material Flow Equipment, LLC sister Company with
Material Flow Solutions, Inc
7010 NW 23rd Way, Ste A
Gainesville, FL 32653 USA
Direct: 352-303-9123

Distributed by:



Particulate Systems
a Division of Micromeritics
4356 Communications Drive
Norcross, GA 30093 USA
US Sales: 770-662-3633
International Sales: 770-662-3600

Fax: 770-662-3696